


UK-437719

Chemical Properties

CAS No. : 856708-54-8
Formula: C₂₉H₄₁F₂N₅O₂
Molecular Weight: 529.677
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	UK-437719 is a bioactive chemical.
-------------	------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8879 mL	9.4397 mL	18.8793 mL
5 mM	0.3776 mL	1.8879 mL	3.7759 mL
10 mM	0.1888 mL	0.944 mL	1.8879 mL
50 mM	0.0378 mL	0.1888 mL	0.3776 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481