


Colistimethate Sodium

Chemical Properties

CAS No. :	30387-39-4
Formula:	C57H108N16O25S4
Molecular Weight:	1545.82
Appearance:	no data available
Storage:	keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Colistin-1,3,5,9-tetramethanesulfonic A acid is used for the treatment of serious infections due to selected aerobic Gram-negative pathogens in patients with limited treatment options.
Targets(IC50)	Others

Solubility Information

Solubility	H2O: 10mg/mL (6.47mM) DMSO: 60 mg/mL (38.81 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6469 mL	3.2345 mL	6.4691 mL
5 mM	0.1294 mL	0.6469 mL	1.2938 mL
10 mM	0.0647 mL	0.3235 mL	0.6469 mL
50 mM	0.0129 mL	0.0647 mL	0.1294 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Wilkinson RE, et al. Outcome analysis of colistin-treated burn center patients. Burns. 2017 Apr 11.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481