


Vedroprevir

Chemical Properties

CAS No. : 1098189-15-1
 Formula: C₄₅H₆₀ClN₇O₉S
 Molecular Weight: 910.53
 Appearance: no data available
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Vedroprevir, also known as GS9451, is a selective hepatitis C virus (HCV) NS3 protease inhibitor in the development of the treatment of genotype 1 (GT1) HCV infection.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0983 mL	5.4913 mL	10.9826 mL
5 mM	0.2197 mL	1.0983 mL	2.1965 mL
10 mM	0.1098 mL	0.5491 mL	1.0983 mL
50 mM	0.022 mL	0.1098 mL	0.2197 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Meissner EG, Kohli A, Virtaneva K, Sturdevant D, Martens C, Porcella SF, McHutchison JG, Masur H, Kottlilil S.
 Achieving sustained virologic response after interferon-free hepatitis C virus treatment correlates with hepatic

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481