


W 6211

Chemical Properties

CAS No. : 22759-46-2
Formula: C₁₄H₂₂N₂O₂
Molecular Weight: 250.34
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	W 6211 is a biochemical.
-------------	--------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9946 mL	19.9728 mL	39.9457 mL
5 mM	0.7989 mL	3.9946 mL	7.9891 mL
10 mM	0.3995 mL	1.9973 mL	3.9946 mL
50 mM	0.0799 mL	0.3995 mL	0.7989 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Moorman JR, Yee R, Bjornsson T, Starmer CF, Grant AO, Strauss HC. pKa does not predict pH potentiation of sodium channel blockade by lidocaine and W6211 in guinea pig ventricular myocardium. J Pharmacol Exp Ther. 1986 Jul;

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481