


W7477

Chemical Properties

CAS No. : 34258-79-2
Formula: C₁₆H₁₄N₂O₄S
Molecular Weight: 330.36
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	W7477 is a biochemical.
-------------	-------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.027 mL	15.135 mL	30.270 mL
5 mM	0.6054 mL	3.027 mL	6.054 mL
10 mM	0.3027 mL	1.5135 mL	3.027 mL
50 mM	0.0605 mL	0.3027 mL	0.6054 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Kinawi A, Siebler D. [The binding of oxacam derivatives to human serum albumin]. Arzneimittelforschung. 1988 Aug;38(8):1089-92. German. Erratum in: Arzneimittelforschung 1988 Dec;38(12):1849. PubMed PMID: 3264167.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481