


YM 534

Chemical Properties

CAS No. : 95923-66-3
Formula: C₂₄H₃₆N₂O₄
Molecular Weight: 416.55
Appearance: no data available
Storage: keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	YM 534 has antitumor activity and platelet aggregation activity, inhibits photoaffinity labeling of deglycosylated P-glycoprotein, and can be used to study leukemia.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4007 mL	12.0034 mL	24.0067 mL
5 mM	0.4801 mL	2.4007 mL	4.8013 mL
10 mM	0.2401 mL	1.2003 mL	2.4007 mL
50 mM	0.048 mL	0.2401 mL	0.4801 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Kanamaru R, et al. Mechanism of action of the synthetic imidazole compound YM 534 on human promyelocytic leukemia line HL-60 cells. Anticancer Drug Des. 1989 Aug;4(2):137-44.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481