


ZM 260384

Chemical Properties

CAS No. : 161229-62-5
Formula: C₁₅H₁₁F₄N₃O₄
Molecular Weight: 373.26
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	ZM 260384 is a potassium channel opener.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6791 mL	13.3955 mL	26.791 mL
5 mM	0.5358 mL	2.6791 mL	5.3582 mL
10 mM	0.2679 mL	1.3395 mL	2.6791 mL
50 mM	0.0536 mL	0.2679 mL	0.5358 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Shaw-Reid CA, Feuston B, Munshi V, Getty K, Krueger J, Hazuda DJ, Parniak MA, Miller MD, Lewis D. Dissecting the effects of DNA polymerase and ribonuclease H inhibitor combinations on HIV-1 reverse-transcriptase activities.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481