


Potassium phosphate dibasic

Chemical Properties

CAS No. :	7758-11-4
Formula:	HK ₂ O ₄ P
Molecular Weight:	174.174
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Potassium phosphate dibasic is a useful organic compound for research related to life sciences. The catalog number is T35382 and the CAS number is 7758-11-4.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.7415 mL	28.7076 mL	57.4152 mL
5 mM	1.1483 mL	5.7415 mL	11.483 mL
10 mM	0.5742 mL	2.8708 mL	5.7415 mL
50 mM	0.1148 mL	0.5742 mL	1.1483 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481