

Data Sheet (Cat.No.T35471)

A 410099.1 amide-PEG3-amine

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Functionalized IAP ligand for PROTAC research and development; incorporates an IAP ligand plus a short amide-PEG linker with terminal amine ready for conjugation to a target protein ligand. Part of a range of functionalized tool molecules for PROTAC R&D. PROTAC is a registered trademark of Arvinas Operations, Inc., and is used under license.
-------------	---

Solubility Information

Solubility	DMSO: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481