


SU5201

Chemical Properties

CAS No. : 114727-43-4
Formula: C₁₅H₉Cl₂NO
Molecular Weight: 290.14
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	SU5201 (NSC 247030) is a biochemical.
Targets(IC50)	Others,Interleukin

Solubility Information

Solubility	DMSO: 50 mg/mL (172.33 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4466 mL	17.2331 mL	34.4661 mL
5 mM	0.6893 mL	3.4466 mL	6.8932 mL
10 mM	0.3447 mL	1.7233 mL	3.4466 mL
50 mM	0.0689 mL	0.3447 mL	0.6893 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481