


dTAGV-1-NEG

Chemical Properties

CAS No. : 2451573-87-6
Formula: C₆₈H₉₀N₆O₁₄S
Molecular Weight: 1247.56
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Negative control for dTAGV-1.
-------------	-------------------------------

Solubility Information

Solubility	DMSO: 124.75 mg/mL(100 mM) (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8016 mL	4.0078 mL	8.0156 mL
5 mM	0.1603 mL	0.8016 mL	1.6031 mL
10 mM	0.0802 mL	0.4008 mL	0.8016 mL
50 mM	0.016 mL	0.0802 mL	0.1603 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481