


PAR 4 (1-6)

Chemical Properties

CAS No. : 225779-44-2
Formula: C₂₈H₄₁N₇O₉
Molecular Weight: 619.676
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	PAR 4 (1-6) can be used for relevant research in the life sciences. Its product number is T36293 and CAS number is 225779-44-2.
-------------	---

Solubility Information

Solubility	DMF: 50 mg/mL DMSO: 50 mg/mL PBS (pH 7.2): 10 mg/mL Ethanol: 50 mg/mL (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6137 mL	8.0687 mL	16.1374 mL
5 mM	0.3227 mL	1.6137 mL	3.2275 mL
10 mM	0.1614 mL	0.8069 mL	1.6137 mL
50 mM	0.0323 mL	0.1614 mL	0.3227 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481