# Data Sheet (Cat.No.T36468)



#### ARN14988

#### **Chemical Properties**

CAS No.: 1502027-70-4

Formula: C16H24ClN3O5

Molecular Weight: 373.83

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

### **Biological Description**

Description ARN14988 is a potent inhibitor of acid ceramidase (IC50 = 12.8 nM for the human enzyme). It inhibits acid ceramidase activity and increases levels of C16 dihydro ceramide and C16 ceramide in A375, G361, M14, MeWo, MNT-1, and SK-MEL-28 melanoma cells. ARN14988 also reduces growth of A375 and G361 melanoma cells (EC50s = 41.8 and 67.7  $\mu$ M, respectively).

## **Solubility Information**

Solubility DMSO: 30 mg/mL DMF: 30 mg/mL

Ethanol: Slightly soluble

DMSO:PBS (pH 7.2) (1:2): 0.33 mg/mL

(< 1 mg/ml refers to the product slightly soluble or insoluble)

#### **Preparing Stock Solutions**

	1mg	5mg	10mg
1 mM	2.675 mL	13.3751 mL	26.7501 mL
5 mM	0.535 mL	2.675 mL	5.350 mL
10 mM	0.2675 mL	1.3375 mL	2.675 mL
50 mM	0.0535 mL	0.2675 mL	0.535 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com