

Cal Green™ 1 (potassium salt)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Appearance: no data available

Storage: keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Cal Green 1 is a cell-impermeant fluorescent calcium indicator that is characterized by high quantum yield and low phototoxicity. Its peak excitation and emission wavelengths (506 and 531 nm, respectively) are comparable to standard fluorescein dyes, making Cal Green 1 appropriate for fluorescent microscopy. Cal Green 1 is ~5-fold brighter than fluo-3 at saturating calcium levels.
-------------	---

Solubility Information

Solubility	H2O: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481