

AL 8810 ethyl amide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	AL 8810 is an 11 β -fluoro analog of prostaglandin F2 α (PGF2 α) which acts as a potent and selective antagonist at the FP receptor. AL 8810 ethyl amide is an analog of AL 8810 in which the C-1 carboxyl group has been modified to an N-ethyl amide. This modification is analogous to the PG N-ethyl amides, as typified by Bimatoprost, that have been introduced as alternative PG ocular hypotensive prodrugs. In contrast to AL 8810 which contracted the cat iris, AL 8810 ethyl amide showed no contraction activity at concentrations up to 10 ⁻⁴ M and did not antagonize the activity of PGF2 α -ethanolamide in this system.
-------------	---

Solubility Information

Solubility	Ethanol:PBS (pH 7.2)(1:2): 0.3 mg/mL Ethanol: 15 mg/mL DMF: 15 mg/mL DMSO: 10 mg/mL (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481