


Darutigenol

Chemical Properties

CAS No. : 5940-00-1
Formula: C₂₀H₃₄O₃
Molecular Weight: 322.48
Appearance: no data available
Storage: store at low temperature, keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Darutigenol has obvious antithrombotic effect, its mechanism may be related to inhibition of platelet aggregation and adhesion.
Targets(IC ₅₀)	Others

Solubility Information

Solubility	DMSO: 60 mg/mL(186.1 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.101 mL	15.5048 mL	31.0097 mL
5 mM	0.6202 mL	3.101 mL	6.2019 mL
10 mM	0.3101 mL	1.5505 mL	3.101 mL
50 mM	0.062 mL	0.3101 mL	0.6202 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Chen FF, et al. Molecularly imprinted polymer for the specific solid-phase extraction of kirenol from Siegesbeckia pubescens herbal extract. Talanta. 2012 Jan 30;89:505-12.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481