


Fmoc-Gly-Ser(psi(Me,Me)pro)-OH

Chemical Properties

CAS No. : 1095952-22-9
Formula: C₂₃H₂₄N₂O₆
Molecular Weight: 424.45
Appearance: no data available
Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Fmoc-Gly-Ser(psi(Me,Me)pro)-OH is a dipeptide.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.356 mL	11.780 mL	23.5599 mL
5 mM	0.4712 mL	2.356 mL	4.712 mL
10 mM	0.2356 mL	1.178 mL	2.356 mL
50 mM	0.0471 mL	0.2356 mL	0.4712 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481