


Thalidomide-O-amido-C3-PEG3-C1-NH2

Chemical Properties

CAS No. :	1799711-29-7
Formula:	C27H35F3N4O11
Molecular Weight:	648.589
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Thalidomide-O-amido-C3-PEG3-C1-NH2 is a synthesized conjugate compound that serves as an E3 ligase ligand-linker. It features a cereblon ligand based on Thalidomide and a 3-unit PEG linker. This compound is specifically designed for use in PROTAC technology.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5418 mL	7.709 mL	15.4181 mL
5 mM	0.3084 mL	1.5418 mL	3.0836 mL
10 mM	0.1542 mL	0.7709 mL	1.5418 mL
50 mM	0.0308 mL	0.1542 mL	0.3084 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Lohbeck J, et al. Practical synthesis of a phthalimide-based Cereblon ligand to enable PROTAC development. *Bioorg Med Chem Lett.* 2016;26(21):5260-5262. doi:10.1016/j.bmcl.2016.09.048

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481