


cis-Isolimonenol

Chemical Properties

CAS No. : 22972-51-6
Formula: C₁₀H₁₆O
Molecular Weight: 152.23
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Cis-Isolimonenol ((1S,4R)-p-Mentha-2,8-dien-1-ol) is a constituent of essential oil.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.569 mL	32.845 mL	65.6901 mL
5 mM	1.3138 mL	6.569 mL	13.138 mL
10 mM	0.6569 mL	3.2845 mL	6.569 mL
50 mM	0.1314 mL	0.6569 mL	1.3138 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Vinutha M, et al. EFFECT OF DIFFERENT ABIOTIC STRESS ON ESSENTIAL OIL YIELD FROM AERIAL PART OF CYMBOPOGON FLEXUOSUS (NEES EX STEUD) WATS. Int J Approx Reason 1.982. 2016 Dec;4(12):1490-1504.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481