


(S)-ErSO

Chemical Properties

CAS No. :	2407860-34-6
Formula:	C22H13F6NO3
Molecular Weight:	453.33
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	(S)-ErSO is an inactive dextrorotatory enantiomer of ErSO.
Targets(IC50)	Others
In vitro	(S)-ErSO is inactive in MCF-7 cells while (R)-ErSO induces dose-dependent death[1].

Solubility Information

Solubility	DMSO: 180 mg/mL (397.1 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2059 mL	11.0295 mL	22.059 mL
5 mM	0.4412 mL	2.2059 mL	4.4118 mL
10 mM	0.2206 mL	1.1029 mL	2.2059 mL
50 mM	0.0441 mL	0.2206 mL	0.4412 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

David J. Shapiro, et al. Activators of the unfolded protein response. WO2020009958A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481