# Data Sheet (Cat.No.T39992)



## FGFR1 inhibitor-2

Chemical Proper	ies
CAS No. :	2410612-08-5
Formula:	C25H22F5N3O3
Molecular Weight:	507.461
Appearance: 🦲	no data available
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year

### **Biological Description**

Description	FGFR1 inhibitor-2 (with an IC50 of 4.55 μM in MDA-MB-231 cells) is a potent inhibitor of FGFR1. It is specifically useful for studying metastatic triple-negative breast cancer.
In vitro	FGFR1 inhibitor-2, when applied at concentrations of 3 and 6 µM to MDA-MB-231 cells for a duration of 48 hours, effectively diminishes FGFR1 expression, as evidenced by Western Blot Analysis[1]. Additionally, a dose-dependent increase in the apoptotic index of MDA-MB-231 cells is observed, reaching a significant 10.6-fold elevation at a concentration of 5 µM, indicating a pronounced pro-apoptotic effect at this dosage[1].

#### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9706 mL	9.853 mL	19.706 mL
5 mM	0.3941 mL	1.9706 mL	3.9412 mL
10 mM	0.1971 mL	0.9853 mL	1.9706 mL
50 mM	0.0394 mL	0.1971 mL	0.3941 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

#### Reference

Ashraf-Uz-Zaman M, et al. Design, synthesis and structure-activity relationship study of novel urea compounds as FGFR1 inhibitors to treat metastatic triple-negative breast cancer. Eur J Med Chem. 2021;209:112866.

#### Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481