


24:0 Lyso PC

Chemical Properties

CAS No. : 325171-59-3
 Formula: C32H66NO7P
 Molecular Weight: 607.854
 Appearance: no data available
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	24:0 Lyso PC, a lysophospholipid (LyP), has potential applications in mRNA drug delivery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6451 mL	8.2257 mL	16.4514 mL
5 mM	0.329 mL	1.6451 mL	3.2903 mL
10 mM	0.1645 mL	0.8226 mL	1.6451 mL
50 mM	0.0329 mL	0.1645 mL	0.329 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481