Data Sheet (Cat.No.T40710)



2-NP-AHD

Chemical Propert	ies	
CAS No. :	623145-57-3	o II
Formula:	C10H8N4O4	HN
Molecular Weight:	248.198	o Nan
Appearance:	no data available	
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year	

Biological Description	C.	
Description	2-NP-AHD is a 2-nitrophenyl derivative of AHD, which is a metabolite of nitrofurans, a	
	class of antibiotics. It serves as an indicator for the detection of illicit use of nitrofuran	
	drugs.	

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.029 mL	20.145 mL	40.2901 mL
5 mM	0.8058 mL	4.029 mL	8.058 mL
10 mM	0.4029 mL	2.0145 mL	4.029 mL
50 mM	0.0806 mL	0.4029 mL	0.8058 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Li Z, et, al. Simultaneous detection of four nitrofuran metabolites in honey by using a visualized microarray screen assay. Food Chem. 2017 Apr 15;221:1813-1821.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481