

Sephadex LH 20

Chemical Properties

CAS No. : 9041-37-6

Formula:

Molecular Weight:

Sephadex LH 20

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Sephadex LH 20 is a suitable medium for the isolation of natural compounds found in various sources, including red wine and pigments, aiding in their purification and extraction processes.
-------------	--

Reference

Cyril Auger, et al. The red wine extract-induced activation of endothelial nitric oxide synthase is mediated by a great variety of polyphenolic compounds. Mol Nutr Food Res. 2010 Jul;54 Suppl 2:S171-83.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481