

Sephadex G 50

Chemical Properties

CAS No. : 9048-71-9

Formula:

Molecular Weight:

Sephadex G 50

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Sephadex G 50 is a gel filtration medium commonly utilized in gel permeation chromatography for effective fractionation of glycopeptide mixtures.
-------------	---

Reference

Galpin IJ, et, al. Gel filtration of protected peptides on sephadex G-50 in hexamethylphosphoramide containing 5% water. J Chromatogr. 1975 Mar 19;106(1):125-9.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481