

Gel filtration medium G-100

Chemical Properties

CAS No. : 9050-94-6

Formula:

Molecular Weight:

Gel filtration medium G-100

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Gel filtration medium G-100, a substance designed for protein purification, enables the separation and purification of proteins through gel filtration.
-------------	---

Reference

Andrews P. Estimation of the molecular weights of proteins by Sephadex gel-filtration. Biochem J. 1964 May;91(2): 222-33.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481