# Data Sheet (Cat.No.T4499)



## Mequitazine

## **Chemical Properties**

CAS No.: 29216-28-2 Formula: C20H22N2S

Molecular Weight: 322.47

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

# **Biological Description**

Description	Mequitazine (Virginan) is an effective, nonsedative and long-acting histamine H1 antagonist.
Targets(IC50)	Histamine Receptor
In vitro	Mequitazine is widely studied and used for allergic disorders such as hay fever and urticaria[1]. Mequitazine demonstrates significant bactericidal effects against all the tested clinical isolates including Ps. aeruginosa. Its effect against the Gram-positive isolates is more pronounced [2].
In vivo	Mequitazine antagonizes the effect of histamine in guinea-pig ileum competitively.  Mequitazine at 107 produces a parallel shift of the dose-response curve to acetylcholine in the rat duodenum. Mequitazine at highest concentration shows anticholinergic activity [3]. Mequitazine inhibits contractile responses to KCl, phenylephrine (PE), 5-hydroxytryptamine (5-HT), and Ca2+ in rat aorta [4].

## **Solubility Information**

		_
Solubility	DMSO: 50 mg/mL (155.05 mM),	
	(< 1 mg/ml refers to the product slightly soluble or insoluble)	

## **Preparing Stock Solutions**

	1mg	5mg	10mg
1 mM	3.1011 mL	15.5053 mL	31.0106 mL
5 mM	0.6202 mL	3.1011 mL	6.2021 mL
10 mM	0.3101 mL	1.5505 mL	3.1011 mL
50 mM	0.062 mL	0.3101 mL	0.6202 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Page 1 of 2 www.targetmol.com

## Reference

Gonnot V, et al. Expedient synthesis of mequitazine an antihistaminic drug by palladium catalyzed allylic alkylation of sodium phenothiazinate. Chem Pharm Bull (Tokyo). 2009 Nov;57(11):1300-2.



Page 2 of 2 www.targetmol.com