


2,3-Diaminopropionic acid hydrochloride

Chemical Properties

CAS No. : 1482-97-9
Formula: C₃H₉ClN₂O₂
Molecular Weight: 140.57
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

HCl



Biological Description

Description	2,3-Diaminopropionic acid hydrochloride (3-Amino-L-alanine hydrochloride) is a competitive cystathionase(CTH) inhibitor.
Targets(IC50)	Others

Solubility Information

Solubility	H ₂ O: 100 mg/mL(711.39 mM),Sonication is recommended. DMSO: insoluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.1139 mL	35.5695 mL	71.1389 mL
5 mM	1.4228 mL	7.1139 mL	14.2278 mL
10 mM	0.7114 mL	3.5569 mL	7.1139 mL
50 mM	0.1423 mL	0.7114 mL	1.4228 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Mushahwar I K , Koeppel R E . Rat liver L-diaminopropionate ammonia lyase. Identification as cystathionase[J]. Journal of Biological Chemistry, 1973, 248(21):7407-7411.
Sondhi S M , Singh N , Johar M , et al. Synthesis,

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481