


Clopirac

Chemical Properties

CAS No. : 42779-82-8
Formula: C₁₄H₁₄ClNO₂
Molecular Weight: 263.72
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Clopirac is a potent and orally active prostaglandin synthetase inhibitor that is an anti-inflammatory agent [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7919 mL	18.9595 mL	37.919 mL
5 mM	0.7584 mL	3.7919 mL	7.5838 mL
10 mM	0.3792 mL	1.896 mL	3.7919 mL
50 mM	0.0758 mL	0.3792 mL	0.7584 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481