Data Sheet (Cat.No.T60647)



A4B17

Chemical Properties

CAS No.:

Formula: C14H7F4NS

Molecular Weight: 297.27

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description A4B17 is an inhibitor of androgen receptor N-terminal with potential for androgen-

responsive prostate cancer treatment.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3639 mL	16.8197 mL	33.6395 mL
5 mM	0.6728 mL	3.3639 mL	6.7279 mL
10 mM	0.3364 mL	1.682 mL	3.3639 mL
50 mM	0.0673 mL	0.3364 mL	0.6728 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com