


LP-284

Chemical Properties

CAS No. : 2412580-47-1
Formula: C₁₆H₂₀N₂O₄
Molecular Weight: 304.34
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	LP-284 kills solid tumours. LP-284 is a potent DNA alkylating agent that can be used in the hematologic cancers studies with compromised DNA repair, for example, mantle cell lymphoma (MCL) [1] [2].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2858 mL	16.429 mL	32.858 mL
5 mM	0.6572 mL	3.2858 mL	6.5716 mL
10 mM	0.3286 mL	1.6429 mL	3.2858 mL
50 mM	0.0657 mL	0.3286 mL	0.6572 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481