


Tropatepine

Chemical Properties

CAS No. :	27574-24-9
Formula:	C ₂₂ H ₂₃ NS
Molecular Weight:	333.49
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Tropatepine is an orally active antagonist of anticholinergic muscarinic. Tropatepine is an antiparkinsonian agent that can be used to combat against extrapyramidal syndrome caused by neuroleptic drugs [1][2].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9986 mL	14.993 mL	29.9859 mL
5 mM	0.5997 mL	2.9986 mL	5.9972 mL
10 mM	0.2999 mL	1.4993 mL	2.9986 mL
50 mM	0.060 mL	0.2999 mL	0.5997 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481