


IU1

Chemical Properties

CAS No. : 314245-33-5
Formula: C₁₈H₂₁FN₂O
Molecular Weight: 300.37
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	IU1 (IC ₅₀ =4.7 μM), a reversible, specific human USP14 proteasome inhibitor, can penetrate the cell.
Targets(IC ₅₀)	DUB, Autophagy
In vitro	IU1 binds specifically to the activated form of USP14. IU1 can potentially inhibit USP14 by preventing its docking on the proteasome, exhibiting little or no activity toward 8 other DUBs, IsoT, UCH37, BAP1, UCH-L1, UCH-L3, USP15, USP2, USP7. USP14 inhibition is rapidly established upon addition of IU1 and rapidly reversed upon its removal. IU1 inhibits USP14 induced chain trimming and decreases electrophoretic mobility of Ub-CCNB species. IU1 enhances proteasomal degradation of Ub-CCNB in the presence of USP14. IU1 promotes degradation of tau and depletes TDP-43, ATXN3, and glial fibrillary acidic protein (GFAP) in proteotoxic mechanisms. [1]
Kinase Assay	High-throughput screening: Screening is conducted at the ICCB-Longwood screening facility. 10 μL of recombinant USP14 protein are dispensed into each well of a 384-well low volume plate in duplicate, using a Wellmate plate dispenser. 33.3 nL of compound from the library are pin-transferred into the wells using a Seiko pin transfer robotic system, followed by pre-incubation for about 30 min. The last two columns of each plate are used for positive and negative controls for the assay. To initiate the enzyme reaction, 10 μL of VS-proteasome plus Ub-AMC mixture are added to each well, using a Wellmate dispenser. Samples are then incubated for another 45 min. Ub-AMC hydrolysis is measured at Ex355/Em460 using an Envision plate reader. The final concentrations of USP14, VS-proteasome and Ub-AMC are 15 nM, 1 nM and 0.8 μM, respectively. The final concentration of test compound is approximately 17 μM. Enzymes and substrates are prepared in Ub-AMC assay buffer (50 mM Tris-HCl (pH 7.5), 1 mM EDTA, 1 mM ATP, 5 mM MgCl ₂ , 1 mM DTT, and 1 mg/mL ovalbumin).
Cell Research	MTT (Only for Reference)

Solubility Information

Solubility	DMSO: 30 mg/mL (100 mM), Ethanol: 15 mg/mL (50 mM),
------------	--

(< 1 mg/ml refers to the product slightly soluble or insoluble)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3292 mL	16.6461 mL	33.2923 mL
5 mM	0.6658 mL	3.3292 mL	6.6585 mL
10 mM	0.3329 mL	1.6646 mL	3.3292 mL
50 mM	0.0666 mL	0.3329 mL	0.6658 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Byung-Hoon L, et al. Nature, 2010, 467(7312), 179-184.

Wu W, Xu H, Liao C, et al. Blockade of USP14 potentiates type I interferon signaling and radiation-induced

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481