


Propiverine

Chemical Properties

CAS No. : 60569-19-9
Formula: C₂₃H₂₉NO₃
Molecular Weight: 367.48
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Propiverine is a potent antimuscarinic agent that effectively inhibits cellular calcium influx, leading to the reduction of muscle spasm. It exhibits neurotropic and musculotropic effects specifically on the smooth muscles of the urinary bladder. Propiverine proves valuable in research concerning overactive bladder (OAB) [1] [2].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7212 mL	13.6062 mL	27.2124 mL
5 mM	0.5442 mL	2.7212 mL	5.4425 mL
10 mM	0.2721 mL	1.3606 mL	2.7212 mL
50 mM	0.0544 mL	0.2721 mL	0.5442 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481