


Antiparasitic agent-5

Chemical Properties

CAS No. :	2494276-55-8
Formula:	C ₂₀ H ₁₆ ClN ₃ O ₃
Molecular Weight:	381.81
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Antiparasitic agent-5 (compound 8h) exhibits potent selectivity as an antiparasitic agent against <i>Leishmania infantum</i> (L. infantum), with an IC ₅₀ value of 2.50 μM. Additionally, this compound displays cytotoxic effects on HepG2 cells with a CC ₅₀ value of 6.78 μM [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6191 mL	13.0955 mL	26.191 mL
5 mM	0.5238 mL	2.6191 mL	5.2382 mL
10 mM	0.2619 mL	1.3096 mL	2.6191 mL
50 mM	0.0524 mL	0.2619 mL	0.5238 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481