


IDO-IN-16

Chemical Properties

CAS No. : 2663591-36-2
Formula: C₂₂H₂₁F₃N₄O
Molecular Weight: 414.42
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	IDO-IN-16 (compound 5) is an IDO inhibitor (IC ₅₀ : 36 nM).
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.413 mL	12.0651 mL	24.1301 mL
5 mM	0.4826 mL	2.413 mL	4.826 mL
10 mM	0.2413 mL	1.2065 mL	2.413 mL
50 mM	0.0483 mL	0.2413 mL	0.4826 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481