


AMPK activator 4

Chemical Properties

CAS No. : 2493239-46-4
Formula: C₂₄H₂₁ClN₂O₃
Molecular Weight: 420.89
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	AMPK activator 4 is a potent and selective AMPK activator that does not inhibit mitochondrial complex I. AMPK activator 4 selectively activates AMPK in muscle tissue and dose-dependently improves glucose tolerance in normal mice, significantly reduces fasting glucose levels and improves insulin resistance in db/db diabetic mice. AMPK activator 4 has hypoglycemic effects. activator 4 has hypoglycemic effects.
Targets(IC50)	AMPK
In vitro	AMPK activator 4 (compound B10) (10-20 µM; 72 h) slightly decreases the viability of HepG2 cells.[1] AMPK activator 4 (0-20 µM; 12-24 h) induced dose-dependent phosphorylation of ACC in C2C12 myotube cells, and the phosphorylation of ACC was highly consistent with the phosphorylation of AMPK. Increased the phosphorylation of AMPK and ACC in human hepatocarcinoma HepG2 and HuH-7 cells.[1]
In vivo	AMPK activator 4 (100 mg/kg; Intragastrical administration; daily for 9 weeks; db/db mice) significantly lowered the fasting blood glucose in db/db mice (2-6 weeks). Significantly improved insulin resistance in db/db mice (for 9 weeks). The body weight and the serum levels of alanine aminotransferase (ALT) and aspartate aminotransferase (AST) (two markers of liver injury) were not significantly affected. It shows an anti-hyperglycemic effect.[1]

Solubility Information

Solubility	DMSO: 50 mg/mL (118.8 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3759 mL	11.8796 mL	23.7592 mL
5 mM	0.4752 mL	2.3759 mL	4.7518 mL
10 mM	0.2376 mL	1.188 mL	2.3759 mL
50 mM	0.0475 mL	0.2376 mL	0.4752 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Sun G, et al. Discovery of AdipoRon analogues as novel AMPK activators without inhibiting mitochondrial complex I. Eur J Med Chem. 2020;200:112466.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481