


Basoparib

Chemical Properties

CAS No. : 1858179-75-5
Formula: C₁₈H₂₁F₂N₇O₃
Molecular Weight: 421.4
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Basoparib is an inhibitor of ribose polymerase (PARP) and has shown antitumour effects.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.373 mL	11.8652 mL	23.7304 mL
5 mM	0.4746 mL	2.373 mL	4.7461 mL
10 mM	0.2373 mL	1.1865 mL	2.373 mL
50 mM	0.0475 mL	0.2373 mL	0.4746 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481