


IDO-IN-14

Chemical Properties

CAS No. : 2568302-02-1
Formula: C₂₃H₂₃ClN₄O₂
Molecular Weight: 422.91
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	IDO-IN-14 is an IDO inhibitor (IC ₅₀ : 0.6928 nM).
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3646 mL	11.8228 mL	23.6457 mL
5 mM	0.4729 mL	2.3646 mL	4.7291 mL
10 mM	0.2365 mL	1.1823 mL	2.3646 mL
50 mM	0.0473 mL	0.2365 mL	0.4729 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481