


Sunobinop

Chemical Properties

CAS No. : 1126793-40-5
Formula: C₂₆H₃₃N₃O₃
Molecular Weight: 435.56
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Sunobinop (S 117957) is an opioid receptor-like orphan receptor (ORL1) modulator.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2959 mL	11.4795 mL	22.9589 mL
5 mM	0.4592 mL	2.2959 mL	4.5918 mL
10 mM	0.2296 mL	1.1479 mL	2.2959 mL
50 mM	0.0459 mL	0.2296 mL	0.4592 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481