


HPK1-IN-11

Chemical Properties

CAS No. : 2734167-68-9
Formula: C₂₇H₂₅N₅O₂
Molecular Weight: 451.52
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	HPK1-IN-11 is a potent HPK1 inhibitor. HPK1-IN-11 has potential for the study of HPK1-related diseases.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2147 mL	11.0737 mL	22.1474 mL
5 mM	0.4429 mL	2.2147 mL	4.4295 mL
10 mM	0.2215 mL	1.1074 mL	2.2147 mL
50 mM	0.0443 mL	0.2215 mL	0.4429 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481