


Erteberel

Chemical Properties

CAS No. :	533884-09-2
Formula:	C ₁₈ H ₁₈ O ₃
Molecular Weight:	282.33
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Erteberel (LY500307) is a potent, selective estrogen receptor β agonist with EC ₅₀ of 0.66 nM, 32-fold selectivity against estrogen receptor α . Phase 2.
Targets(IC ₅₀)	Estrogen Receptor/ERR, Estrogen/progestogen Receptor

Solubility Information

Solubility	Ethanol: 57 mg/mL (201.89 mM), DMSO: 57 mg/mL (201.89 mM), H ₂ O: <1 mg/mL, < 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.542 mL	17.7098 mL	35.4195 mL
5 mM	0.7084 mL	3.542 mL	7.0839 mL
10 mM	0.3542 mL	1.771 mL	3.542 mL
50 mM	0.0708 mL	0.3542 mL	0.7084 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Norman BH, et al. J Med Chem, 2006, 49(21), 6155-6157.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481