

ATM Inhibitor-2

Chemical Properties

CAS No. :

Formula: C₂₆H₃₁N₇O₃

Molecular Weight: 489.57

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	ATM Inhibitor -2 is a potent and selective inhibitor of ATM (IC ₅₀ <1 nM).
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0426 mL	10.213 mL	20.4261 mL
5 mM	0.4085 mL	2.0426 mL	4.0852 mL
10 mM	0.2043 mL	1.0213 mL	2.0426 mL
50 mM	0.0409 mL	0.2043 mL	0.4085 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 36 Washington Street, Wellesley Hills, MA 02481