


PD-1/PD-L1-IN-25

Chemical Properties

CAS No. : 2768759-52-8
Formula: C₂₆H₂₄ClN₃O₅
Molecular Weight: 493.94
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	PD-1/PD-L1-IN-25 is an inhibitor of PD-1/PD-L1 interaction (IC ₅₀ : 16.17 nM) that effectively activates anti-tumor immunity of T cells in PBMCs. PD-1/PD-L1-IN-25 can be used to study cancer.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0245 mL	10.1227 mL	20.2454 mL
5 mM	0.4049 mL	2.0245 mL	4.0491 mL
10 mM	0.2025 mL	1.0123 mL	2.0245 mL
50 mM	0.0405 mL	0.2025 mL	0.4049 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481