


NSC 23766 trihydrochloride

Chemical Properties

CAS No. :	1177865-17-6
Formula:	C ₂₄ H ₃₅ N ₇ ·3HCl
Molecular Weight:	530.96
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	NSC 23766 trihydrochloride (Rac1 Inhibitor) is an inhibitor of Rac GTPase targeting Rac activation by GEFs (IC ₅₀ : ~50 μM); no inhibitory for the closely related targets, RhoA or Cdc42.
Targets(IC ₅₀)	Apoptosis,Rho,Ras
In vitro	NSC23766 is identified to fit into a surface groove of Rac1 known to be critical for GEF specification. NSC23766 effectively inhibits Rac1 binding and activation by the Rac-specific GEF Trio or Tiam1 in a dose-dependent manner without interfering with the closely related Cdc42 or RhoA binding or activation by their respective GEFs or with Rac1 interaction with BcrGAP or effector PAK1. [1] NSC 23766 is active in regulating Rac GTPase functions on cytoskeleton and many cell functions including cell cycle, cell growth, adhesion, migration and gene transcription. NSC 23766 (50 μM) potentially blocks serum or platelet-derived growth factor-induced Rac1 activation and lamellipodia formation without affecting the activity of endogenous Cdc42 or RhoA in NIH 3T3 cells. NSC 23766 reduces Trio or Tiam1 but not Vav, Lbc, Intersectin, or a constitutively active Rac1 mutant-stimulated NIH 3T3 cells growth and suppresses Trio, Tiam1, or Ras-induced cell transformation. NSC23766 dose-dependently inhibits PC-3 cells proliferation and anchorage-independent growth. 25 μM NSC23766 inhibits the PC-3 cell invasion through Matrigel by 85%. [1] 50 μM NSC 23766 inhibits thrombin-induced activation of Rac1 and Rac2 in human platelets, as well as platelet aggregation. [2] NSC23766 prevents Aβ ₄₀ and Aβ ₄₂ production in swAPP-HEK293cells without affecting Notch and sAPPα. NSC23766 prevents γ-secretase activity in cell, but not act as a direct γ-secretase inhibitor. NSC23766 dose-dependently reduces levels of secreted and intracellular Aβ ₄₀ with IC ₅₀ of 48.94 μM. 50 μM NSC 23766 inhibits release of Aβ ₄₂ by 57.97%. [3]
In vivo	NSC23766 induces mobilization of hematopoietic stem cells/progenitors. Intraperitoneal administration of NSC23766 (2.5 mg/kg) into the "poorly mobilizing" C57Bl/6 mouse strain leads to a two-fold increase in circulating hematopoietic stem cells/progenitors 6 hr after injection. [2] NSC23766 alleviates lipopolysaccharide-induced acute pulmonary injury in mice. Treatment with NSC23766 at 1 or 3 mg/kg not only reduces the inflammatory cells infiltration and MPO activities, but also inhibits pro-inflammatory mediators, tumor necrosis factor-α and interleukin-1β, mRNA expression. NSC23766 also reduces Evans Blue and albumin accumulation in LPS-challenged lungs. [6]

A DRUG SCREENING EXPERT

Kinase Assay	Rho GTPase activity assay: Cells are grown in log phase in a 10-cm dish, and are starved in 0.5% serum medium or indicated otherwise for 24 h before lysis in a buffer containing 20 mM Tris HCl (pH 7.6), 100 mM NaCl, 10 mM MgCl ₂ , 1% Nonidet P-40, 10% glycerol, and 1× protease inhibitor mixture. Lysates are clarified, the protein concentrations are normalized, and the GTP-bound Rac1 in the lysates is measured by an effector domain pull-down assay. For the His6-PAK1 PBD pull-down assay, cell lysates are incubated with Ni ²⁺ -agarose-immobilized His6-PAK1 PBD domain (~ 1 µg each) purified from E. coli for 30 min. The Ni ²⁺ -agarose co-precipitates are washed twice in the wash buffer and analyzed by immunoblotting with anti-Rac1 monoclonal antibody.
Cell Research	Cells (1.5 × 10 ⁴ /mL) are seeded in each well of 96-well tissue culture plates with 200 µL of medium. After 24 hours of plating, the medium is replaced with 200 µL of fresh medium containing NSC23766 at the indicated concentrations. At the end of the treatment period 20 µL of MTS solution are added to each well and incubated at 37 °C for 2 hours. Absorbance at 490 nm is read on a 96-well plate reader.(Only for Reference)

Solubility Information

Solubility	H ₂ O: 53.1 mg/mL (100 mM), DMSO: 53.1 mg/mL (100 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8834 mL	9.4169 mL	18.8338 mL
5 mM	0.3767 mL	1.8834 mL	3.7668 mL
10 mM	0.1883 mL	0.9417 mL	1.8834 mL
50 mM	0.0377 mL	0.1883 mL	0.3767 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Ma X, Tan X, Yu B, et al. DOCK2 regulates antifungal immunity by regulating RAC GTPase activity. Cellular & Molecular Immunology. 2022: 1-17.
Gao Y, et al. Proc Natl Acad Sci U S A, 2004, 101(20), 7618-7623.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481