


STING agonist-10

Chemical Properties

CAS No. : 2259308-69-3
Formula: C₂₅H₂₀ClF₄N₃O₂
Molecular Weight: 505.89
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	STING agonist-10 is a potent activator of the STING small molecule cyclic urea class (EC ₅₀ : 2600 nM).STING activation is a highly promising immunotherapy.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9767 mL	9.8836 mL	19.7671 mL
5 mM	0.3953 mL	1.9767 mL	3.9534 mL
10 mM	0.1977 mL	0.9884 mL	1.9767 mL
50 mM	0.0395 mL	0.1977 mL	0.3953 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481