


HDAC/HSP90-IN-3

Chemical Properties

CAS No. : 2700035-54-5
Formula: C₂₆H₃₃N₅O₆
Molecular Weight: 511.57
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	HDAC/HSP90-IN-3 is a selective and potent dual inhibitor of fungal Hsp90 (IC ₅₀ : 0.83 μ M) and HDAC (IC ₅₀ : 0.91 μ M), exhibiting antifungal effects against azole-resistant <i>Candida albicans</i> . HDAC/HSP90-IN-3 exhibits inhibitory effects on important virulence factors and is able to downregulate the resistance genes ERG11 and CDR1.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9548 mL	9.7738 mL	19.5477 mL
5 mM	0.391 mL	1.9548 mL	3.9095 mL
10 mM	0.1955 mL	0.9774 mL	1.9548 mL
50 mM	0.0391 mL	0.1955 mL	0.391 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481