


CX-6258 hydrochloride hydrate

Chemical Properties

CAS No. : 1353858-99-7
Formula: C₂₆H₂₇Cl₂N₃O₄
Molecular Weight: 516.42
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	CX-6258 hydrochloride hydrate is a potent and selective inhibitor of pan-Pim kinase that acts on Pim-1 (IC ₅₀ : 5 nM), Pim-2 (IC ₅₀ : 25 nM) and Pim-3 (IC ₅₀ : 16 nM).
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9364 mL	9.682 mL	19.3641 mL
5 mM	0.3873 mL	1.9364 mL	3.8728 mL
10 mM	0.1936 mL	0.9682 mL	1.9364 mL
50 mM	0.0387 mL	0.1936 mL	0.3873 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481