Data Sheet (Cat.No.T63678)



RORyt modulator 2

Chemical Propert	ties	
CAS No. :	2247083-47-0	H,C M
Formula:	C28H29ClN2O4S	
Molecular Weight:	525.06	
Appearance: 🦲	no data available	
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year	

Biological Description	Le la
Description	RORyt modulator 2 is a modulator of retinol-related orphan receptor yt (RORyt) (IC50<50 nM) and can be used in the study of RORyt-mediated diseases such as pain, inflammation, asthma, rheumatoid arthritis, COPD, multiple sclerosis, psoriasis, colitis, neurodegenerative diseases and cancer.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9045 mL	9.5227 mL	19.0454 mL
5 mM 📀	0.3809 mL	1.9045 mL	3.8091 mL
10 mM	0.1905 mL	0.9523 mL	1.9045 mL
50 mM	0.0381 mL	0.1905 mL	0.3809 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481