


AEBSF hydrochloride

Chemical Properties

CAS No. :	30827-99-7
Formula:	C ₈ H ₁₁ ClFNO ₂ S
Molecular Weight:	239.69
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	AEBSF hydrochloride (Pefabloc SC) is an irreversible inhibitor of serine proteases, such as chymotrypsin, kallikrein, thrombin, plasmin, and trypsin.
Targets(IC ₅₀)	Thrombin,Serine Protease,Influenza Virus
In vitro	In a mouse model of cockroach allergen-induced asthma, AEBSF was able to mitigate respiratory reactions and potential inflammation. Intraperitoneal injection of 76.8 mg/kg AEBSF daily in mice lethally inoculated with Toxoplasma gondii extended their survival period.
In vivo	AEBSF is a serine protease inhibitor that suppresses the lysis of leukemia cells by macrophages without inhibiting the secretion of TNF- α and IL-1 β by these cells. In five different human cell lines, AEBSF inhibits the formation of amyloid-beta (A β) by blocking the secretion of both neural and non-neural β -secretase products. Additionally, AEBSF disrupts embryonic growth within the endometrium and modifies protein secretion mechanisms to prevent the adhesion of HeLa cells in human umbilical vein endothelial cells.
Cell Research	The HeLa cells suspended in RPMI-1640 media containing 10% FCS are plated into each well of a 96-well microplate (5 $\times$ 10 ³ cells/200 μ L/well). After incubation for 24 h at 37°C, cells are treated with different doses of AEBSF (0, 25, 50, 100 μ g/mL) for 48 h. Then, 20 μ L fresh 3-(4,5)-dimethylthiazoliazol-2-yl-3,5-diphenyltetrazolium bromide (MTT) reagent (5 μ g/ μ L) is added into each well, and cells are cultured at 37°C in 5% CO ₂ for another 4 h. The media are discarded carefully, and 150 μ L DMSO is added. Absorbance is read on a microplate reader at dual wavelengths of 540 and 620 nm.

Solubility Information

Solubility	H ₂ O: 24 mg/mL (100 mM), DMSO: 60 mg/mL (250.32 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1721 mL	20.8603 mL	41.7206 mL
5 mM	0.8344 mL	4.1721 mL	8.3441 mL
10 mM	0.4172 mL	2.086 mL	4.1721 mL
50 mM	0.0834 mL	0.4172 mL	0.8344 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Citron M, et al. Neuron. 1996, 17(1), 171-179.

Jiang T Y, Pan Y F, Wan Z H, et al. PTEN status determines chemosensitivity to proteasome inhibition in